

Using Relative Humidity (RH) Sensing Technology; A Solvent-Free Alternative To Karl Fischer Titration



The Health Care industry has increased its needs for specialized devices over the past decade, which has led to a new frontier of resin and polymer development designed to keep the quality of care high while minimizing cost.

With these goals in mind, the resins being used for medical devices are scrutinized more thoroughly than other resins that require less regulatory compliance. Analysing a product for outgassing, deformation, and reactivity, among other things, has become part of the daily routine for manufacturers, moulders, and final inspection personnel before an item can be shipped or used. This additional testing and control also includes the amount of water that is allowed in the resins, since this will greatly influence the final product's rigidity, consistency, and lifetime, as well as the quality of care that will be provided to the customer. Oftentimes the quality control of the materials is closely monitored using testing equipment defined in an IQ/OQ/PQ: installation qualification (IQ), operational qualification (OQ), and performance qualification (PQ) to ensure that the instruments are effective, and the quality of the product is consistent.

Moisture Determination

As an alternative to the Coulometric Karl Fisher titration, Relative Humidity (RH) sensor moisture detection was first used as a method for determination of water in materials in 1997, with the introduction of the Computrac® 3000 Moisture Analyser by Arizona Instrument LLC. This method uses a thermoset polymer capacitor that has a selective response when in the presence of water, the same way that many RH sensors work in traditional settings such as houses, laboratory controlled environments, and dry boxes. Medical device resins are sealed in a sample vial, and then transported into an oven chamber with inert gas blown through it. As the material gets hot, water molecules evolve off and are carried to the sensor via the carrier gas. The sensor is exposed to the water molecules and a measurable change in the electronic activity takes place. This method requires no solvents, making it an environmentally friendly

Registered Address

ABLE Instruments & Controls Ltd
Cutbush Park, Danehill, Lower Earley,
Reading, Berkshire, RG6 4UT. UK.

Reading

Tel: +44 (0)118 9311188 | Email: info@able.co.uk

Aberdeen

Tel: +44 (0)1224 725999 | Email: ab@able.co.uk

Web

able.co.uk

E-commerce

247able.com



alternative to traditional chemical titration. The instrument provides in-situ moisture measurements, which allows users to monitor performance in real time. Additionally, it has detection Limit of 10ppm, and is more rugged than Karl Fisher titrators, making it a suitable instrument for moisture analysis in manufacturing facilities, as well as Quality Control and inspection labs.

This technology is now being adopted as the standard test method and is described by ASTM D7191, Standard Test Method for Determination of Moisture in Plastics by Relative Humidity Sensor. This instrument also meets the high demands of performance given in an IQ/OQ/PQ.

With major advances in technology, the medical device community is also taking advantage of new RAPID loss-on-drying methods for moisture determination. These instruments use the same principle as traditional loss-on-drying techniques, but address the shortcomings of the method.

Sample material is heated on a balance and real time measurements are providing immediate feedback and moisture concentration. The Computrac[®] MAX[®] 4000XL instrument, manufactured by Arizona Instrument LLC, provides a parameter development expert program that allows users to optimize testing conditions, such as sample size, test ending criteria, testing temperature, idle temperature, temperature rate, etc. The chassis of this instrument is made of steel, which prevents cracking in the case and cool air from entering the testing chamber, which would influence the results.

These new techniques are being adopted as standard testing methods and are described by ASTM D6980 -12, Standard Test Method for Determination of Moisture in Plastics by Loss in Weight. Like the Vapor Pro[®] 3100L, the MAX[®] 4000XL meets the performance standards set forth in typical IQ/OQ/PQ testing.

Testing

Sample Prep – Medical grade thermoplastic polyurethane (TPU), polycarbonate (PC), and nylon 6/6 resins were selected for analysis. The materials were stored wet in a 1

gallon plastic Ziploc bags prior to testing. An initial analysis of TPU was conducted to determine the water content prior to drying. The material was then dried in the Dri-Air HP4-X 25 plastics drying hopper for 6 hours prior to testing. The material remained in the dryer during testing due to the hygroscopic properties of the material. Once it was determined that 6 hours was sufficient for testing, the PC and nylon 6/6 were also dried for 6 hours. These materials were stored in Mason jars, upside down, to prevent head-space moisture from influencing the results.

Test Conditions – Reference testing was conducted using the Mitsubishi CA-100 Coulometric Karl Fischer titrator.

The parameters for TPU testing were:

sample size – 0.5g +/- 0.1g,
temperature – 90°C,
purge/preheat/cooling – 1/2/2,
ending sensitivity – 0.1µg/sec.
PC: sample size – 0.3g +/- 0.1g,
temperature – 220°C,
purge/preheat/cooling – 1/2/2,
ending sensitivity – 0.1µg/sec.
Nylon6/6: sample size – 0.5g +/- 0.1g,
temperature – 220°C,
purge/preheat/cooling – 1/2/2,
ending sensitivity – 0.1µg/sec.

Corollary testing was conducted using the Computrac[®] Vapor Pro[®] 3100L.

The parameters for TPU testing were:

sample size – 2g +/- 0.2g,
temperature – 105°C,
purge – 50 seconds,
ending criteria – rate<0.1µg/sec.
PC: sample size – 1.0g +/- 0.1g,
temperature – 180°C,
purge – 40 seconds,
ending criteria – rate<0.30µg/sec.
Nylon 6/6: sample size – 1.0g +/- 0.1g,
temperature – 240°C, purge – 30 seconds,
ending criteria – rate<0.05µg/sec.

Registered Address

ABLE Instruments & Controls Ltd
Cutbush Park, Danehill, Lower Earley,
Reading, Berkshire, RG6 4UT. UK.

Reading

Tel: +44 (0)118 9311188 | Email: info@able.co.uk

Aberdeen

Tel: +44 (0)1224 725999 | Email: ab@able.co.uk

Web

able.co.uk

E-commerce

247able.com



Results

	Karl Fischer		Vapor® Pro 3100L	
	Result (ppm)	Test time	Result (ppm)	Test time
	47	6:53	51	12:32
	58	7:30	74	14:36
	67	7:04	53	12:20
	79	7:46	78	15:08
			69	13:18
Average	62.8		65.0	
S.D.	13.6		11.0	
RSD	21.6		16.9	

Table 1. Comparative results of TPU testing


	Karl Fischer		Vapor® Pro 3100L	
	Result (ppm)	Test time	Result (ppm)	Test time
	48	2:00	53	3:20
	51	2:00	58	4:34
	52	2:00	60	5:22
	54	2:00	54	4:24
	46	2:00	69	5:44
Average	50		59	
S.D.	3.0		6.0	
RSD	6.0		10.2	

Table 2. Comparative results for Lexan 1 testing

	Karl Fischer		Vapor® Pro 3100L	
	Result (ppm)	Test time	Result (ppm)	Test time
	352	6:53	444	10:54
	419	7:30	410	10:40
	502	7:04	442	11:28
	486	7:46	477	10:24
	483		504	11:12
Average	448		455	
S.D.	63.0		36	
RSD	14.1		7.91	


Table 3. Comparative results for Polyone testing

Wet TPU Total Moisture Curve




Graph 1. Total moisture curve of pre-paid TPU

Dry TPU Total Moisture Curve



Graph 2. Total moisture curve of TPU dried for 6 hours at 200°C

Rate Graph for Dried PC Resin



Graph 3. Real time rate of loss for dried PC material

Registered Address

ABLE Instruments & Controls Ltd
Cutbush Park, Danehill, Lower Earley,
Reading, Berkshire, RG6 4UT. UK.

Reading

Tel: +44 (0)118 9311188 | Email: info@able.co.uk

Aberdeen

Tel: +44 (0)1224 725999 | Email: ab@able.co.uk

Web

able.co.uk

E-commerce

247able.com



From the tables, the results using the two different instruments with similar testing conditions correlate to each other for all three medical grade materials. For two of the three materials the Vapor Pro® did show a statistical improvement in the relative standard deviation, but did require a slightly longer test time than the Karl Fischer. Additionally, the Vapor Pro® provided real time data points that could be used to graph the total moisture and rate change curves. This allows for better monitoring of the product, or diagnosing possible problems with the instrument, and allows users to adjust testing parameters. Examining the rate curve for PC, it was noticed that the rate had a significant drop at $0.30\mu\text{g}\cdot\text{sec}^{-1}$. This is what was used to determine the ending criteria for this material. This feature was not available for Karl Fischer titrator.

Conclusion

The development of an alternative to Karl Fischer moisture analyser has been achieved, and can be used for moisture specific analysis of medical device grade resins. The Computrac® Vapor Pro® 3100L moisture analyser successfully uses Relative Humidity sensor technology for selective and accurate moisture measurement. The instrument reduces the use of hazardous organic solvents makes it an environmentally friendly alternative, when compared to current Karl Fischer technology. The results between the two methods of detection of H₂O content for three separate resins, TPU, PC, and nylon 6/6 strongly correlate, with the Vapor Pro® 3100L providing real-time data that can be used to provide a complete moisture profile of the materials.

References available on request



Registered Address

ABLE Instruments & Controls Ltd
Cutbush Park, Danehill, Lower Earley,
Reading, Berkshire, RG6 4UT. UK.

Reading

Tel: +44 (0)118 9311188 | Email: info@able.co.uk

Aberdeen

Tel: +44 (0)1224 725999 | Email: ab@able.co.uk

Web

able.co.uk

E-commerce

247able.com

